

EXAM ROOM SUPPLIES

# KINDERGARTEN HANGING EYE CHART

### Features

- Shows easy to distinguish objects
- 20' distance
- Eyelet allows chart to be hung on wall
- This is a chart for measuring visual acuity at a distance - in other words, how well you can see things that are far away

	overall Width		overall Height		test Distance		UOM
	inch	cm	inch	cm	ft	m	
1243	11	27.9	22	55.9	20	6.1	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2005, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a registered trademark of GF Health Products, Inc.