

## STANDARD WHEELCHAIRS

# TRAVELER® L3 PLUS WHEELCHAIR

### L3 PLUS SERIES FIXED BACK WHEELCHAIRS

Everest & Jennings® introduces the Traveler® L3 Plus as the lightweight wheelchair addition to the Traveler® line, a brand respected for performance and durability. The Traveler® L3 Plus lightweight frame weighs less than 36 pounds making it the ideal K0003 Medicare coded wheelchair built for short and long term use as well as rental needs.

This lightweight wheelchair comes standard with non-marking polyurethane tires, industry compatible front rigging and adjustable anti-tippers along with a detachable 2" seat extension with hardware. It features a dual position axle, multi-position caster forks for 18 or 20 inch, 16.3" or 17 1/4" seat height adjustment and offers flip-back desk, adjustable height and full arm options for easy lateral movement. The Traveler® L3 Plus is manufactured with the same commitment to quality Everest & Jennings® is known for.

#### Features

- Weighs less than 36 lb without front rigging
- Available in 16", 18", 20" seat widths including 2" removable seat extension kit converting the chair from 16" depth to 18" depth if necessary
- Two position axle and caster fork provides 18" to 20" seat height in regular models and seat height 16.3" to 17.9" in the XVI models
- Hammertone finish and carbon steel frame provide strength and durability
- Black nylon upholstery with standard chart pocket on back
- Includes fastener tab under seat connecting both seat and back upholstery to eliminate gap
- Anti-tippers come standard
- Non-marking polyurethane tires with 24" rear wheels with 8" casters in the regular models and 20" rear wheels with 6" casters in the XVI models
- Tool-free adjustable front rigging with heel loops standard on swingaway footrests
- Flip-back padded arms available in fixed desk length, adjustable height desk and detachable flip-back full length
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0003

#### Warranty:

Five (5) year limited warranty



*fixed-back model with flip-back desk arms (shown) and adjustable swingaway footrest (shown) — other configurations available*





**Traveler® L3 Plus Series Fixed-Back Wheelchair**

Item	HCPCS+ code	Overall Chair Width		Overall Chair Depth		Overall Chair Height		Folded Chair Width		Seat Width		Seat Depth		Seat Height		Back Height		Wheel Diameter		Caster Diameter		Product Weight <sup>8</sup> w/o front rigging		Weight <sup>1</sup> Capacity EVENLY DISTRIBUTED		Flip-Back Full Arms	Flip-Back Desk	Flip-Back Adjust. Height Desk	Detach. Flip-Back Full Fixed Desk	Flip-Back Adjust. Height Desk	Swingaway Footrest	Detach. Swingaway Footrest	Elevating Legrest	XVI Model	Rear Anti-Tippers		
		inch	cm <sup>1</sup>	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kgf	lb	kgf											ams	padded arms
3F012120	+E0951	26.00	66.04	29.00	73.66	35.50	90.17	12.25	31.12	18.00	45.72	16.00 <sup>3</sup>	40.64 <sup>3</sup>	18	45.72	17.25	43.82	24.00	60.96	8.00	20.32	36	16.3	300	136.08	—	—	—	—	—	—	—	—	—	—	—	—
3F012130	+E0990													18	45.72											—	—	—	—	—	—	—	—	—	—	—	
3F012140	+E0951													18	45.72											—	—	—	—	—	—	—	—	—	—	—	—
3F012220	+E0951	24.00	60.96	29.00	73.66	35.50	90.17	12.25	31.12	16.00	40.64	16.00 <sup>3</sup>	40.64 <sup>3</sup>	18	45.72	17.25	43.82	24.00	60.96	8.00	20.32	36	16.3	300	136.08	—	—	—	—	—	—	—	—	—	—	—	—
3F012230	+E0990													18	45.72											—	—	—	—	—	—	—	—	—	—	—	—
3F012320	+E0951+E2201	28.00	71.12	29.00	73.66	35.50	90.17	12.25	31.12	20.00	50.80	16.00 <sup>3</sup>	40.64 <sup>3</sup>	18	45.72	17.25	43.82	24.00	60.96	8.00	20.32	36	16.3	300	136.08	—	—	—	—	—	—	—	—	—	—	—	—
3F012330	+E2201+E0990													18	45.72											—	—	—	—	—	—	—	—	—	—	—	—
3F012340	+E0951+E2201													18	45.72											—	—	—	—	—	—	—	—	—	—	—	—
3F013120	+E0951+E0973	26.00	66.04	29.00	73.66	35.50	90.17	12.25	31.12	18.00	45.72	16.00 <sup>3</sup>	40.64 <sup>3</sup>	18	45.72	17.25	43.82	24.00	60.96	8.00	20.32	36	16.3	300	136.08	—	—	—	—	—	—	—	—	—	—	—	
3F013220	+E0951+E0973+E2201	24.00	60.96							16.00	40.64			18	45.72											—	—	—	—	—	—	—	—	—	—	—	
3F013320	+E0951+E0973+E2201	28.00	71.12							20.00	50.80			16.30	41.40											—	—	—	—	—	—	—	—	—	—	—	
3F014120	+E0951	25.78	65.48	20.00	50.80	33.00	83.82	12.25	31.12	18.00	45.72	16.00 <sup>3</sup>	40.64 <sup>3</sup>	16.30	41.40	17.25	43.82	20.00	50.80	6.00	15.24	31	14.1	300	136.08	—	—	—	—	—	—	—	—	—	—	—	
3F014220	+E0951	23.81	60.48							16.00	40.64			17.90	45.47											—	—	—	—	—	—	—	—	—	—	—	
3F014320	+E0951+E2201	27.75	70.49							20.00	50.80															—	—	—	—	—	—	—	—	—	—	—	

■ Standard ■ Optional ■ Available ■ Accessory — Not Apply / Available

<sup>1</sup> designed to meet the requirements of <sup>†</sup> standard rounding method used <sup>1</sup> maximum weight capacity - evenly distributed - refer to user manual <sup>3</sup> +2 inches (5.08 cm) with seat depth extension bar (included) <sup>8</sup> approximate measurement: product weight — as seat width and rigging change, so does product weight ★ options available at additional cost - accessories sold separately

**Parts & Accessories:** ★ sold separately

	HCPCS+	UOM	Item No.		HCPCS+	UOM	Item No.		HCPCS+	UOM	Item No.
Footrest Heel Loop	E0951	pair	90763064(A)	Flip-Back Full Arm Assembly	K0015 +K0019	pair	90763232A	Combination O2 & IV Holder	E2208	each	B10320F
Wheel Lock Fixed Arm (with rubber tips)	—	pair	90763070A	Arm Pad for Full Arm - Pad Only	K0019	pair	90763023A	O2 Cylinder Holder	E2208	each	B10325
Positioning Belt (auto-style, push button) (adjustable length: 27 - 47 inches (68.58 - 119.38 cm))	E0978	each	300-505	2" Seat Extension for 16" seat width	—	kit	90763SE16	O2 Holder Bag (nylon)	E2208	each	EJ02BAG07
Adjust. Height Flip-Back Desk Arm Assembly	E0973	pair	90763227	2" Seat Extension for 18" seat width	—	kit	90763SE18	Adjustable Anti-Tipper with Wheel	E0971	pair	90763230A
Flip-Back Desk Arm Assembly	—	pair	90763231A	2" Seat Extension for 20" seat width	—	kit	90763SE20	Adjustable Anti-Tipper with Wheel for XVI	E0971	pair	90763233A
				Rail Plug (dust cap) for use with Seat Extension	—	2 pair	907633CSE	Wheelchair Positioning Tray	—	each	EJTRAY18



GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2019, GF Health Products, Inc. All Rights Reserved. Everest & Jennings is a registered trademark of GF Health Products, Inc.



GF1600013 Rev H124