

TRAPEZE / FLOOR STANDS

PORTABLE PATIENT HELPER TRAPEZE

The Gendron[®] Portable Patient Helper Trapeze unit is an important patient room accessory designed to help the users change positions while in bed and to aid in transferring from bed to chair with easy attendant assistance.

Features

- Durable, robust steel frame
- Helps to position patient in bed and transfer from bed to other equipment with minimal assistance
- Handle bar location fully adjustable using strap with buckle to maximize patient accessibility
- Comes with 2" diameter x 1" wide casters and floor stand base that folds for transportability
- Limited 2 Year Limited Warranty



Specifications:

	overall Height		overall Width (frame and leg width)		overall Depth (leg length)		Boom Length (extends horizontally over bed)		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
GH935850	72	182.88	36	91.44	63	160.02	35	88.90	154	69.85	850	385.55	1/ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.gendroninc.com | www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved.
Gendron is a registered trademark of GF Health Products, Inc.