


## DIAGNOSTIC INSTRUMENTS

# ORTHOPEDIC GONIOMETER

### Features

- The orthopedic goniometer is used to measure a joint's range of motion (ROM). There are two "arms"—one that is stationary and one that is movable—that are hinged together
- Plastic



	Length		UOM
	inch	cm	
12-1000	12	30.5	1 ea
12-1001	8	20.3	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a trademark of GF Health Products, Inc.