

DIAGNOSTIC INSTRUMENTS

# STANDARD OTOSCOPE SET

Graham-Field's Standard Otoscope Set offers an LCD screen indicating brightness level and battery life, a durable, ergonomic handle that is easy to control in your hand and comfortable to use, and a hard storage case for storing or transporting.

### Features

- LCD screen indicating the brightness level and battery
- Durable ergonomic handle
- LED bulb with dimmable rheostat for brightness control
- Fiber Optic Illumination
- 3x Magnification
- 25 Lumen
- 2 AA Batteries (Not Included)
- One (1) Year Limited Warranty



		shipping Weight		UOM
		lb	Kg	
1227-1	Standard Otoscope Set	1	0.45	1 set



GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2019, GF Health Products, Inc. All Rights Reserved. Graham-Field is a trademark of GF Health Products, Inc.