

DIAGNOSTIC INSTRUMENTS

# STANDARD OPHTHALMOSCOPE SET

Graham-Field's Standard Ophthalmoscope Set offers an LCD screen indicating brightness level and battery life, a durable, ergonomic handle that is easy to control in your hand and comfortable to use, and a hard storage case for storing or transporting

### Features

- LCD screen indicating the brightness level and battery
- LED bulb with dimmable rheostat for brightness control
- Range of Dioptic Lenses
- 5 different apertures
- 14 Lumen
- 2 AA Batteries (Not Included)
- One (1) Year Limited Warranty



		shipping Weight		UOM
		lb	Kg	
1227-2	Standard Ophthalmoscope Set	1	0.45	1 set



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2019, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a trademark of GF Health Products, Inc.