

DIAGNOSTIC INSTRUMENTS

COMPLETE DIAGNOSTIC KIT

Graham-Field's Complete Diagnostic Kit comes complete for evaluation of eyes, nose, and throat. The durable, ergonomic handle equipped with an LCD screen indicating brightness level and battery life is designed to fit all 6 devices. The LED Bulb's brightness is controlled with a dimmable rheostat.

Included

- Otoscope head with 4 specula: 2.4mm, 3mm, 4mm, 5mm
- Ophthalmoscope head with 5 different apertures
- Metal Nasal Speculum 0-25 mm
- Laryngeal Mirror with 20.33 mm replacement mirror
- Tongue Depressor Holder with plastic depressors
- Operating Otoscope with rotatable lens and 5 optional specula



		shipping Weight		UOM
		lb	Kg	
1227-4	Complete Diagnostic Set	2	0.91	1 set



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2019, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.