

DIAGNOSTIC INSTRUMENTS

BUCK NEUROLOGICAL HAMMER

Features

- Chrome-plated octagonal handle
- Compartment for stainless steel needle and brush
- Dual Black rubber bumper



		Length		UOM
		inch	cm	
1304		7	17.8	1 ea
1304-1	Telescoping handle	4.5 - 7	11.4 - 17.8	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.