


EXAM ROOM SUPPLIES

SNELLEN TYPE HANGING EYE CHART

Features

- With eyelet to be hung easily
- The Snellen eye chart is the familiar chart with the big letter or symbol at the top and progressively smaller letters at the bottom of the chart down to 20/20 or better
- This is a chart for measuring visual acuity at a distance - in other words, how well you can see things that are far away
- Non- reflective, matte finish with green and red color bars
- Screening device
- 20' distance

| | overall Width | | overall Height | | test Distance | | UOM |
|------|---------------|------|----------------|------|---------------|-----|------|
| | inch | cm | inch | cm | ft | m | |
| 1240 | 11 | 27.9 | 22 | 55.9 | 20 | 6.1 | 1 ea |



GF Health Products, Inc. is an intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2005, GF Health Products, Inc. All Rights Reserved. Graham-Field is a trademark of GF Health Products, Inc.