



EXAM ROOM SUPPLIES

ILLITERATE PLASTIC EYE CHART

Features

- With green and red color bars
- Ideal for use with patients who are illiterate, not familiar with the Roman alphabet, or cannot speak
- This is a chart for measuring visual acuity at a distance - in other words, how well you can see things that are far away
- Eight lines or the letter "E" pointing in different directions, with green and red bars
- Durable, plastic



| | | overall Width | | overall Height | | test Distance | | equivalent Distance | | UOM |
|------|---|---------------|------|----------------|------|---------------|-----|---------------------|-----|------|
| | | inch | cm | inch | cm | ft | m | ft | m | |
| 1241 | Illiterate/Tumbling E Hanging Eye Chart | 11 | 27.9 | 22 | 55.9 | 20 | 6.1 | 20 | 6.1 | 1 ea |
| 1262 | Illuminated Illiterate/Tumbling E Eye Chart | 9 | 22.9 | 14 | 35.6 | 10 | 3 | 20 | 6.1 | 1 ea |



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2005, GF Health Products, Inc. All Rights Reserved. Graham-Field is a trademark of GF Health Products, Inc.