

EXAM ROOM SUPPLIES

FEMALE URINAL

Features

- Durable HDPE (High Density PolyEthylene) in blue color
- Single patient use
- Easy-to-clean
- Withstand temperatures up to 230°F
- Molded handle
- 8½"H x 3½"W x 5½"L
- Graduated in oz/cc up to 28oz / 800ml capacity



| | overall Height | | overall Width | | overall Length | | graduated Volume | | UOM |
|------|----------------|------|---------------|-----|----------------|----|------------------|-----|------|
| | inch | cm | inch | cm | inch | cm | oz | ml | |
| 3214 | 8.5 | 21.6 | 3.5 | 8.9 | 5.5 | 14 | 28 | 800 | 1 ea |



GF Health Products, Inc. is an Intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved. Graham-Field is a trademark of GF Health Products, Inc.