

## EXAM ROOM SUPPLIES

# STAINLESS STEEL BED PAN

### Features

- Seamless, sanitary, highly-polished 18-8 stainless steel ware
- Liquid capacity is 4 quarts
- Reusable
- Autoclavable



		overall Length		overall Width		front Height		rear Height		Volume		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	qrt	L	lb	Kg	
3227	Adult	14	35.6	11.375	28.9	4.375	11.1	2.625	6.7	4	3.8	500	226.8	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

**+1 770.368.4700**

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2010, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a trademark of GF Health Products, Inc.