

EXAM ROOM SUPPLIES

STAINLESS STEEL
MALE URINAL

Features

- Stainless Steel
- Autoclavable
- 9 1/2" x 4 1/2"
- Capacity: 1 1/4 – 1 3/4 quarts



	Height		Diameter		Volume		UOM
	inch	cm	inch	cm	oz	ml	
3230	9.5	24.1	4.5	11.4	1.25 – 1.75	37 – 51.8	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.