

EXAM ROOM SUPPLIES

PLASTIC SUNDRY JARS

Labeled and Unlabeled Plastic Sundry Jars are excellent for storing and dispensing a variety of supplies, making them ideal for hospitals and clinics.

Features

- Available in transparent or opaque high impact polystyrene construction
- Each jar (labeled set) permanently imprinted with blue labels: Applicators, Gauze, Tongue Depressors, Cotton and Bandages
- Each jar includes removable overlapping plastic cover
- Set includes five (5) plastic jars complete with plastic lids
- Reusable
- Not Autoclavable
- Latex free
- One (1) year limited warranty



		Height		Diameter		case Weight		UOM
		inch	cm	inch	cm	lb	kg	
3452	Clear, Labeled	7	17.8	4.25	10.8	2.7	1.22	5 ea/st
3453	Clear, Unlabeled	7	17.8	4.25	10.8	2.7	1.22	5 ea/st
3455	Opaque, Labeled	7	17.8	4.25	10.8	2.7	1.22	5 ea/st



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.