

DIAGNOSTIC INSTRUMENTS

ADJUSTABLE FREQUENCY TUNING FORK

Features

- Fork has adjustable masses that can be moved up and down to change the pitch
- All aluminum tuning fork, with adjustable weight scale in seven frequencies: C128, D144, E160, F170.7, G192, A213.3 and B240 (weights removed)
- Scale letters are stamped on stems



| | tine Width | | overall Length | | UOM |
|------|------------|------|----------------|------|------|
| | inch | cm | inch | cm | |
| 1325 | 0.375 | 0.95 | 8 | 20.3 | 1 ea |



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.