

LUMEX® Specifications and Dimensions

A Graham-Field Brand

ITEM	Overall Height	Overall Width	Overall Depth	Width Between Arms	Arm Height From Floor	Arm Height From Seat	Seat Height	Seat Width	Seat Depth	Depth Fully Reclined In Trendelenburg	Width With Both Side Tables in Upright Position	Shipping Weight	Max. Weight Capacity
330	44"	31"	38"	23"	19 1/2" - 25 1/2"	5"	19"	21 1/2"	20"	-	-	109 lb	250 lb
561	47"	27"	41"	20"	26"	9"	21"	19"	19"	-	-	125 lb	250 lb
565G	48"	28"	36"	19"	28"	8 1/2"	20"	19"	19 1/2"	-	-	80 lb	250 lb
565PG	46"	26"	34"	17 1/2"	26 1/2"	8"	19"	17"	18"	-	-	74 lb	250 lb
565WG	48"	32"	37"	24"	28"	8 1/2"	20"	24"	19 1/2"	-	-	91 lb	400 lb
565DG	48"	28"	36"	19 1/2"	23 1/2" - 29 1/2"	0" - 8 1/2"	22"	19"	19 1/2"	-	-	85 lb	250 lb
565TG	48"	28"	36"	19 1/2"	23 1/2" - 29 1/2"	0" - 8 1/2"	22"	19"	19 1/2"	-	-	90 lb	250 lb
574G	47"	25"	41"	20"	26"	9"	21"	19"	19"	-	-	110 lb	250 lb
574RG	47"	29"	41"	20"	27 1/2"	9"	21"	19"	19"	-	-	110 lb	290 lb
590G	39 1/2"	27 1/2"	31 1/2"	18 1/2"	25 1/2" - 29 1/2"	8 1/2"	17" - 19"	18"	16" - 18"	-	-	52 lb	225 lb
592G	39 1/2"	27 1/2"	31 1/2"	18 1/2"	25" - 33 1/2"	8" - 12"	17" - 19"	18"	16" - 18"	-	-	59 1/2 lb	225 lb
593G	52"	27 1/2"	33"	18 1/2"	25" - 33 1/2"	8" - 12"	17" - 19"	18"	16" - 18"	-	-	68 lb	225 lb
GF4401-1	42"	23 3/4"	26"	18"	31"	6"	26"	18"	18"	-	-	36 lb	250 lb
GF4401-2	42"	23 3/4"	26"	18"	31"	6"	26"	18"	18"	-	-	36 lb	250 lb
566G	50"	31"	43"	21"	27 1/2"	8"	22"	20 1/2"	20"	75"	51"	145 lb	350 lb
566DG	51 1/2"	33"	38"	23"	28"	8"	22"	20 1/2"	19"	76"	53"	135 lb	350 lb
577RG	47"	31"	35"	22"	25 1/2"	8"	22"	21 1/2"	21"	74"	51"	127 lb	290 lb
587W	50"	33"	36"	23"	26 1/2"	8"	21 1/2"	23 1/2"	20 1/2"	-	-	138 1/2 lb	450 lb

Standard Colors



ROSEWOOD (863)



JADE (857)



BLUE RIDGE (427)



DOE SKIN (851)



GRAY TAUPE (858)



ROYAL BLUE (454)



DOVE GRAY (413)

- Upholstery color may vary slightly from brochure.
- Lumex Ortho-Biotic Seating with an FR prefix designates that the seating meets CA133 Fire Code specifications.
- Please contact Graham-Field Customer Service at 800-347-5678 for details on specialty colors and COM orders.

Graham-Field Health Products is not responsible for typographical errors. Packaging, warranties and product are subject to change without notice.