

OVERBED TABLES

ECONOMY OVERBED TABLE, NON-TILT

Features

- Easy to assemble
- Twin-wheel swivel casters provide smooth-rolling action. Two locking and two non-locking casters
- Walnut, wood-grain laminated particle-board core with plastic edge molding
- Chrome plated or Silvervein hammertone finish "H" base and column for aesthetics and durability
- Designed to meet the requirements of HCPCS code: E0274
- Warranty: One (1) Year



GF8902 - Chrome Plated



GF8902S - Silvervein Hammertone



Specifications:

		table Height		table Width		table Depth		base Width		base Depth		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
GF8902	Chrome Plated	28 - 44"	71.12 - 111.76	30"	76.2	15"	38.1	26.2"	66.55	15.25"	38.74	22	9.98	30	13.61	1/ea
GF8902S	Silvervein Hammertone	29.5" - 42.5"	74.93 - 107.95	28"	71.12	15.4"	39.12	26.4"	67.06	15.07"	38.28	20.5	9.3	50	22.68	1/ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2010, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF-100003Rev/D22