



PATIENT LIFTS AND SLINGS

MAXICARE PATIENT LIFT

PLx00 SERIES BARIATRIC PATIENT LIFTING SYSTEM

When transferring a bariatric patient, safety is key for both the patient and caregiver. The Gendron Patient Lift is designed to offer a safer and robust solution to patient transfers.

Features

- Heavy-gauge steel construction
- Four-Point spreader bar with rotation tilts to enhance comfort and safety
- The hand pendant controls all functions of the lift including extending and retracting the legs and the opening and closing of the base
- Lifting and lowering along with the left and right rotational movement are also hand pendant operated
- Lower base height maximizes stability and fits under beds
- 24V DC motor allows quiet smooth operation
- Easy access emergency stop button
- Emergency Power Off button available if needed in event of a power loss
- Powder-coated finish
- Optional digital scale available
- Available in 3 options to support 700 lb (317.5 kg), 800 lb (362.8 kg) and 1,000 (453.5 kg) Weight Capacities

Warranty

- Frame and spreader bar five (5) year limited
- Actuator one (1) year limited

Product Info



Specifications

Item	UOM	Overall Width		Overall Length		Base Width <i>extension legs open</i>				Spreader Bar Height				Floor Clearance		Product ¹ Weight		Weight Capacity EVENLY DISTRIBUTED		Actuator	Frame & Bar
		inch	cm	inch	cm	minimum	maximum	minimum	maximum	minimum	maximum	inch	cm	lb	kg	lb	kg	warranty			
PL700	1/ea	85.00	215.90	61.69	156.69	53	134.62	85.00	215.9	40.81	103.66	81.94	208.13	4.5	11.43	351	159.21	700	317.51	12	60
PL800	1/ea	85.00	215.90	61.69	156.69	53	134.62	85.00	215.9	43.31	110.01	78.56	199.54	4.75	12.07	354	160.57	800	362.87	12	60
PL1000	1/ea	85.00	215.90	61.69	156.69	53	134.62	85.00	215.9	43.31	110.01	78.56	199.54	4.75	12.07	354	160.57	1000	453.59	12	60



[™] GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.gendroninc.com | www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved.
Gendron is a registered trademark of GF Health Products, Inc.