

REPLACEMENT TOPS




REPLACEMENT TOPS

Make your exam tables look new at a fraction of the cost. All components within in the replacement top are all new including wood, foam, vinyl and hardware.

- Industry standard colors *plus*
- Replacement Tops fit leading brand tables



GF100RT-xx
shown in violet | VI
100



GF104RT-xx
shown in teal | TE
104



GF204RT-xx
shown in british tan | BT
204



GF304RT-xx
shown in ocean gray | OG
304



GF404RT-xx
shown in doe | DO
404



GF405RT-xx
shown in shadow | SH
405



GF420RT-xx
shown in artichoke | AR
420



GF460RT-xx
shown in british tan | BT
460

Specifications:

	100 104		204 ¹		304		404 405		420 460	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
Overall Width:	27.50	69.85	27.50	69.85	27.50	69.85	27.50	69.85	27.50	69.85
Overall Length:	55.00	139.70	55.00	139.70	55.00	139.70	60.00	152.40	60.00	152.40
Footrest Width:	20.00	50.80	19.00	48.26	20.00	50.80	20.50	52.07	19.00	48.26
Footrest Length:	18.00	45.72	18.00	45.72	18.00	45.72	17.50	44.45	18.00	45.72
Shipping Weight:	lb	kg	lb	kg	lb	kg	lb	kg	lb	kg
	55	24.95	55	24.95	55	24.95	35	15.88	55	24.95

Features:	100 104	204 ¹	304	404 405	420 460
Hook & Loop Attachment:	—	—	—	included	—
Paper Roll Holder:	included	—	—	—	included
Mounting Hardware:	included	—	—	—	included
Lift Handles:	—	not included ¹	—	—	—
Top Replacement Knobs:	—	not included ¹	—	—	—
Individual Product Links:	🔗 🔗	🔗	🔗	🔗 🔗	🔗 🔗

¹ 2009 year model and later do not include lift handles and top attachment knobs
If the OEM model has a seamless top, request paper roll holder when ordering
— Does Not Apply 🔗 Active Hyperlink

Standard Color Options (-xx)

Available in 24 color choices to complement most decors.



Swatches are digital reproductions. Perceived texture and coloring may vary on finished product. Monitor and/or printer settings will affect appearance of color. Swatch Card is available from your representative.

GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



© 2018, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.