

## STANDARD WHEELCHAIRS

# TRAVELER® SE PLUS WHEELCHAIR

### SE PLUS SERIES FIXED BACK WHEELCHAIRS

The Traveler® SE Plus continues to incorporate the legendary Everest & Jennings® wheelchair frame with saddle-joint tubular construction with an improved weight capacity of 300 lb evenly distributed.

It provides an effective combination of a standard wheelchair for transportation and a two-inch seat height adjustment for a hemi position.

#### Features

- Available in 16", 18" and 20" seat widths, with 16" seat depth
- Black vinyl upholstery with standard chart pocket on back
- Industry compatible tool-free front rigging with standard heel loops on swingaway footrests with composite footplates
- Carbon steel frame with hammertone finish provides strength and durability
- Comfortable desk and full length arms with padded armrests
- Two-position axle and caster fork provide two-inch seat height adjustment for a hemi position
- Maximum Weight Capacity 300 lb -EVENLY DISTRIBUTED
- Meets Medicare codes K0001/K0002



*fixed-back model with removable desk arms (shown) and adjustable swingaway footrest (shown) — other configurations available*

# Traveler® SE Plus Series Fixed-Back Wheelchair

Item	HCPCS <sup>+</sup> code	Overall Chair Width		Overall Chair Depth		Overall Chair Height		Folded Chair Width		Seat Width		Seat Depth		Seat Height		Back Height		Wheel Diameter		Caster Diameter		Product <sup>®</sup> Weight w/o front rigging		Weight <sup>†</sup> Capacity EVENLY DISTRIBUTED		Detach. Full	Detach. Desk	Fixed Full	Swingaway Footrest	Elevating Legrest	Rear Anti-Tippers	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg†	lb	kg†	arms	front rigging							
3E050100	+E0951	26.50	67.31	31.50	80.01	35.50	90.17	12.25	31.12	18.00	45.72	16.00	40.64	18	45.72	15.50	39.37	24.00	60.96	8.00	20.32	36.5	28.58	300	136.08	—	—	■	■	—	■	
3E050110	+K0195													18	45.72												—	—	■	—	■	■
3E050120	+E0951																										—	■	—	■	—	■
3E050130	+K0195																										—	■	—	—	■	■
3E050140	+E0951																										■	—	—	■	—	■
3E050150	+K0195																										■	—	—	—	■	■
3E050220	+E0951	24.50	62.23	31.50	80.01	35.50	90.17	12.25	31.12	16.00	40.64	16.00	40.64	18	45.72	15.50	39.37	24.00	60.96	8.00	20.32	36.5	28.58	300	136.08	—	■	—	■	—	■	
3E050230	+K0195																										—	■	—	—	■	■
3E050240	+E0951																										■	—	—	■	—	■
3E050320	+E0951+E2201	28.50	72.39	31.50	80.01	35.50	90.17	12.25	31.12	20.00	50.80	16.00	40.64	18	45.72	15.50	39.37	24.00	60.96	8.00	20.32	36.5	28.58	300	136.08	—	■	—	■	—	■	
3E050330	+E0951+E2201																										—	■	—	—	■	■
3E050340	+E0951+E2201																										■	—	—	■	—	■
3E050350	+E0951+E2201																										■	—	—	■	—	■

■ Standard ■ Optional ■ Available ■ Accessory — Not Apply / Available

† designed to meet the requirements of † standard rounding method used † maximum weight capacity - evenly distributed - refer to user manual  
 \* approximate measurement: product weight — as seat width and rigging change, so does product weight \* options available at additional cost - accessories sold separately

Product Info

**Parts & Accessories:** \* sold separately

	HCPCS <sup>+</sup>	UOM	Item No.		HCPCS <sup>+</sup>	UOM	Item No.		HCPCS <sup>+</sup>	UOM	Item No.
Footrest Heel Loop	E0951	pair	90763064(A)	Combination O2 & IV Holder	E2208	each	B10320F	Adjustable Anti-Tipper with Wheel	E0971	pair	90763230A
Wheel Lock Lever Extension	E0961	pair	90763062	O2 Cylinder Holder	E2208	each	B10325	Anti-Fold Bar Lower Frame	—	each	9252380
Positioning Belt (auto-style, push button) adjustable length: 27 - 47 inches (68.58 - 119.38 cm)	E0978	each	300-505	O2 Holder Bag (nylon)	E2208	each	EJ02BAG07	Wheelchair Positioning Tray	—	each	EJTRAY18



GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

**+1 770.368.4700**

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2019, GF Health Products, Inc. All Rights Reserved. Everest & Jennings is a registered trademark of GF Health Products, Inc.



GF1900082RevC24