

DIAGNOSTIC INSTRUMENTS

CHROME PLATED EAR SYRINGE

Features

- Brass construction
- Chrome plated
- Includes shield and two tips
- Available in 3 sizes
- Latex free
- Professional use only

One (1) Year Limited Warranty: Product is warranted to be free from defects in material and workmanship. Replacement will be issued for any product which is found defective under the terms of this warranty.



	Volume		Unit				Bulbous Tip				Tapered Tip				UOM
			inner diameter		length		inner diameter		length		inner diameter		length		
	oz	ml	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	
2807	2	59.1	1.122	2.85	6.625	16.83	0.089	0.23	1.906	4.84	0.094	0.24	1.998	5.07	1 ea
2808	3	88.7	1.403	3.56	6.75	17.15	0.094	0.24	1.877	4.77	0.096	0.24	1.951	4.96	1 ea
2809	4	118.3	1.403	3.56	7.75	19.69	0.094	0.24	1.87	4.75	0.0975	0.25	2.003	5.09	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.