

EMERGENCY/FIRST AID

BERMAN-TYPE ORAL AIRWAYS

Berman-Type Oral Airways are made of smooth, flexible white polyethylene and utilize a dual channel design that allows a suction catheter or endotracheal tube to slide into the pharyngeal space even if the airway is partially obstructed. The disposable airways are available in six sizes to provide proper sizing for each patient.

Features

- Reinforced bite block prevents airway collapse
- Disposable for single-patient use
- Non-sterile, non-toxic polyethylene
- Gas/cold sterilizable
- Not made with natural rubber latex
- Made in the U.S.A.



		Length	UOM
		mm	
3786-12#0	Neonatal	40	12/pk
3786-12#1	Infant	50	12/pk
3786-12#3	Small	80	12/pk

Note: Airway will stay in place without securement by keeping patient's head in a tilted position. Taping an oropharyngeal airway in place is contraindicated as it can lead to aspiration if the patient starts to gag or cough and can't expel the device. Airways should be removed intermittently according to your facility's policy and the patient's condition to inspect oral cavities and provide mouth care. Proper sizing is essential to ensure airway patency; device intended to be used by properly trained healthcare professionals.



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2008, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.