

EMERGENCY/FIRST AID

# DISPOSABLE TOURNIQUET

### Features

- Commonly used amongst Phlebotomists they aid the process of obtaining blood samples by restricting the flow of blood temporarily and making veins more visible
- Designed to ensure no any cross-spread between patients

### Latex

- Smooth latex, disposable, easy to apply and release

### Latex-free

- Ideal for use on patients with sensitive skin
- Eliminates risk of allergic reaction to latex protein



		Width		Length		UOM
		inch	cm	inch	cm	
4109	Latex	1	2.5	18	45.7	100/bx
4109-1	Latex	1	2.5	18	45.7	250/bx
4109-1LF	Latex-free	1	2.5	18	45.7	250/bx



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a trademark of GF Health Products, Inc.