

EXAM ROOM SUPPLIES

FOX ALUMINUM EYE SHIELD

Features

- Lightweight 1060 aluminum
- Perforate convex eye shield with slots
- Individually wrapped
- Weight: 5 gm
- Overall dimensions: 7.3 cm x 6 cm
- Thickness: 0.1 cm
- Autoclavable
- Does not contain latex
- Optional elastic banded protective cloth cover



| | | overall Dimensions | Thickness | Weight | UOM |
|-----------|--|-----------------------|-----------|--------|-------|
| | | cm | cm | g | |
| 1276 | Fox Aluminum Eye Shield | 7.3 cm x 6 cm | 0.1 | 5 | 12/bx |
| 1276GA LB | Fox Aluminum Eye Shield With Light Blue Protective Cloth Cover | 7.3 cm x 6 cm | 0.1 | 5 | 50/bx |
| 1276GA W | Fox Aluminum Eye Shield With White Protective Cloth Cover | 7.3 cm x 6 cm | 0.1 | 5 | 50/bx |
| 1276G LB | Replacement Cloth Cover for Fox Aluminum Eye Shield - Light Blue | 7.3 cm x 6 cm | 0.1 | 5 | 50/bx |



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.