

## EXAM ROOM SUPPLIES

# TRAUMA SHEARS

### Features

- 7 ½" all purpose shears
- Easily cut leather, plastic casts and bandages
- Stainless steel with serrated non-slip cutting edge
- Sterilizable and autoclavable
- Designed purposely for external use only



		Length		UOM
		in	cm	
1358-1	Economy, black	7.5	19.1	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

**+1 770.368.4700**

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2008, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a trademark of GF Health Products, Inc.