

EXAM ROOM SUPPLIES

# STAINLESS STEEL UTILITY CART

### Features

- Provides storage space
- Ergonomic handles

### Components

- **Pipe Frame:** Stainless Steel
- **Steel Shelves:** 3 Stainless Steel Shelves
- **Nuts/ Bolts:** Stainless Steel
- **Wheels:** Thermo Plastic Rubber
- **Wheel Cover:** Co-Polymer
- **Caster:** M.S. Chrome Plated
- **Dimensions:** 16"W x 24"L x 30"H (Floor to Top Shelf Height)
- **Overall Dimensions:** 16.25"W x 28.5"L x 36"H
- **Load Capacity:** 400 lb - EVENLY DISTRIBUTED



	overall Width		overall Depth		overall Height		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	lb	Kg	
8146	16.25	41.3	28.5	72.4	36	91.4	400	181.44	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2008, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a trademark of GF Health Products, Inc.