

## LABORATORY

# TEST TUBE BRUSHES

### Features

- Firm white nylon bristles with radial tufted end on sturdy, twisted, galvanized wire handle
- Sized to clean test tubes with a minimum of effort

Brushes are not autoclavable



|      | overall<br>Length |      | bristle<br>Width |     | bristle<br>Length |     | UOM   |
|------|-------------------|------|------------------|-----|-------------------|-----|-------|
|      | inch              | cm   | inch             | cm  | inch              | cm  |       |
| 3403 | 8                 | 20.3 | 0.5              | 1.3 | 3                 | 7.6 | 12 ea |
| 3404 | 9                 | 22.9 | 0.75             | 1.9 | 3.5               | 8.9 | 12 ea |



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a trademark of GF Health Products, Inc.