


LABORATORY

SMALL BORE TEST TUBE BRUSH

Features

- White bristles with ½" sponge end
- Galvanized twisted wire handle
- Ideal for cleaning small culture tubes, test tubes, sahli tubes and all small bore glassware

Brushes are not autoclavable



| | overall Length | | bristle Width | | bristle Length | | UOM |
|------|-------------------|------|------------------|-----|-------------------|-----|-------|
| | inch | cm | inch | cm | inch | cm | |
| 3414 | 6 | 15.2 | 0.5 | 1.3 | 3 | 7.6 | 12 ea |



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.