

LABORATORY

PIPET BRUSH

Features

- Soft white nylon bristles on galvanized wire handle

Brushes are not autoclavable



	overall Length		Large Bristles				Small Bristles				UOM
			Width		Length		Width		Length		
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	
3419	17	43.2	0.75	1.9	3.25	8.3	0.25	0.6	4	10.2	12 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.