

LABORATORY


BIBULOUS PAPER

Features

- Highly absorbent type of paper
- Used for drying slides and cover glasses, absorbing stains
- (50) 4" x 6" (10.2 x 15.2 cm) sheets per book

Technical specifications:

- A-grade cotton fiber
- 0.18mm thick
- 80gm/square meter



	overall Width		overall Length		UOM
	inch	cm	inch	cm	
3718	4	10.2	6	15.2	6 books/pk



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2019, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.