

MEDICAL/SURGICAL | PATIENT ROOM

PATIENT ID BRACELET

Features

- White opaque textured vinyl bracelet with crystal clear card window welded to band plus insert cards
- Adjusts to fit wrist sizes 4" - 9" (10.2 - 22.9 cm)



| | | Wrist Size | | UOM |
|------|-------|------------|-------------|-----------|
| | | inch | cm | |
| 1503 | Adult | 4 - 9 | 10.2 - 22.9 | 500 ea/bx |



-GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2005, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.