


MEDICAL/SURGICAL | PATIENT ROOM

FOLDING PRIVACY SCREENS

Features

- 6 mil white vinyl panels are flame-retardant (NFPA 701 Test Method 1-2010 Edition)
- Easy to set up and to dismantle for storage
- Vinyl panel is washable and replaceable on snap-out rods
- Frame is 1" (25 mm) anodized aluminum tubing
- 2" (51 mm) spherical casters allow for easy movement

To clean vinyl panels, use mild liquid soap and water, rinse.



		UOM
4299C	Replacement (4) Tips and (4) Casters	1 set

	panel Width		panel Height		overall Height		product Weight		UOM
	inch	cm	inch	cm	inch	cm	lb	kg	
4296W	27.5	69.9	56	142.2	70.8	179.8	14.9	6.75	1 ea



GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website.

www.grahamfield.com



© 2011, GF Health Products, Inc. All Rights Reserved. Graham-Field is a trademark of GF Health Products, Inc.