

## SCALES

# PHYSICIAN MECHANICAL BEAM SCALE WITH WHEELS

### Features

- Dual measurement in lb and kg
- Weight capacity up to 490 lb (220 kg), using counterpoise weight of 50 lb (20 kg)
- Weight can be read from the front or back
- Retractable aluminum height rod, Range of 24" – 84" (60 cm – 210 cm)
- Heavy-duty roller wheels allow easy portability of the scale
- Durable enameled steel body
- Cast aluminum beam
- Removable slip-resistant platform cover
- No parts contain latex: Latex Free
- Meeting the OIML, International Organization of Legal Metrology International Recommendation R76
- Cubic Dimensions:  
59.5" x 20.7" x 6.3" (151 cm x 52.5 cm x 16 cm)

**Caution:** Placing the scale on a soft surface such as carpet is not recommended. This could result in inaccurate weight measurement. Ensure scale is on a level, stable and flat surface when in use.

### Warranty

Two (2) year limited warranty on manufacturer defects.

This warranty does not include calibration of the scale. See User Instructions for the details.



	overall Width		overall Depth		overall Height		platform Depth		platform Length		Height Rod (range)		max Weight Capacity		net Weight		gross Weight		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg	lb	kg	lb	kg	
HT485	10.8	27.5	20.9	53	58.5	148.5	14.8	37.5	10.8	27.5	24–84 x 0.25	60–210 x 0.5	490 x 0.25	220 x 0.1	35.5	16	45.2	20.5	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2011, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a trademark of GF Health Products, Inc.