

SPHYGMOMANOMETERS

# AUTOMATIC BLOOD PRESSURE MONITOR

Blood pressure management is one of the most important things you can do to maintain your health and is essential to preventing heart disease and strokes. More doctors are recommending the use of home blood pressure monitors so that patients can regularly, conveniently check their blood pressure.

Graham-Field's Automatic Upper Arm Blood Pressure Monitor makes home blood pressure management stress free with its ease of use and portable design.

### Features

- Automatic Inflation and Deflation Blood Pressure Monitor
- Digital LCD
- Memory feature - 2 user with 90 sets of memories each
- Average value of the latest 3 readings
- Uses 4 AA batteries (not included)
- Includes D-ring cuff designed to be self-applied
- Irregular heartbeat detection



### Specifications:

	overall Width		overall Length		overall Height		UOM
	inch	mm	inch	mm	inch	mm	
1133	4.72	120	4.13	105	2.17	55	1/ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2010, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a trademark of GF Health Products, Inc.