

## SPHYGMOMANOMETERS

# ACCUMAX<sup>TM</sup> COMPLETE INFLATION SYSTEMS

### Features

- Calibrated nylon Accumax<sup>TM</sup> inflation cuff
- Two-tube bladder
- Deluxe bulb & valve



| Item  | Size        | Tubes  | Bladder   | Color | UOM  |
|-------|-------------|--------|-----------|-------|------|
| 2361A | Adult       | 2-Tube | Latex     | Blue  | 1 ea |
| 2361N | Newborn     | 2-Tube | Latex     | Blue  | 1 ea |
| 2361C | Child       | 2-Tube | Latex     | Blue  | 1 ea |
| 2361X | Large Adult | 2-Tube | Latex     | Blue  | 1 ea |
| 2361T | Thigh       | 2-Tube | Latex     | Blue  | 1 ea |
| 2365A | Adult       | 2-Tube | Latex     | Black | 1 ea |
| 2365X | Large Adult | 2-Tube | Latex     | Black | 1 ea |
| 2365T | Thigh       | 2-Tube | Latex     | Black | 1 ea |
| 2360A | Adult       | 2-Tube | Non-Latex | Black | 1 ea |
| 2360C | Child       | 2-Tube | Non-Latex | Black | 1 ea |
| 2360T | Thigh       | 2-Tube | Non-Latex | Black | 1 ea |
| 2360X | Large Adult | 2-Tube | Non-Latex | Black | 1 ea |
| 2366A | Adult       | 2-Tube | Non-Latex | Grey  | 1 ea |
| 2367A | Adult       | 2-Tube | Non-Latex | Blue  | 1 ea |
| 2367C | Child       | 2-Tube | Non-Latex | Blue  | 1 ea |

|          | Limb Range  |             | cuff Width |      | cuff Height |      |
|----------|-------------|-------------|------------|------|-------------|------|
|          | inch        | cm          | inch       | cm   | inch        | cm   |
| Newborn  | 3.75 – 5.25 | 9.5 – 13.3  | 7          | 17.8 | 2           | 5.1  |
| Infant   | 5 – 7.5     | 12.7 – 19.1 | 10.5       | 26.7 | 3           | 7.6  |
| Child    | 7 – 10.5    | 17.8 – 26.7 | 13.5       | 34.3 | 4           | 10.2 |
| Adult    | 9.5 – 15    | 24.1 – 38.1 | 20.5       | 52.1 | 5.5         | 14   |
| LG Adult | 13 – 19     | 33 – 48.3   | 25         | 63.5 | 6.5         | 16.5 |
| Thigh    | 17 – 26     | 43.2 – 66   | 29.5       | 74.9 | 8.5         | 21.6 |



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2022, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a trademark of GF Health Products, Inc.