

## STETHOSCOPES

# LIGHTWEIGHT PANASCOPE<sup>TM</sup> STETHOSCOPES

### Features

- This lightweight aluminum stethoscope offers a dual head chestpiece with a chrome-plated rotating stem
- The non-chill ring is color coordinated with 22" (56.5 cm) latex-free tubing
- Adjustable binaurals are chrome-plated brass with an internal spring
- Comes complete with additional diaphragm and eartips



Item	Size	Color	UOM
500	Adult	Black	1 ea
500GY	Adult	Grey	1 ea
510	Pediatric	Black	1 ea
510GY	Pediatric	Grey	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

**+1 770.368.4700**

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2008, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a trademark of GF Health Products, Inc.