

## STETHOSCOPES

# DELUXE SPRAGUE-RAPPAPORT TYPE PROFESSIONAL STETHOSCOPE

### Features

- Five interchangeable black chestpieces (two diaphragms, three bells) facilitates the detection of high, medium, and low frequency sounds
- Includes three black eartip sets to maximize comfort
- Heavy-walled 22" (56.5 cm) tubing decreases extraneous sounds
- Ideal for infant, pediatric, and adult use
- Black metal binaural
- Accessory pouch included
- Latex free

Item	Color	UOM
602BK	Midnight Black	1 ea



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.grahamfield.com](http://www.grahamfield.com)



© 2009, GF Health Products, Inc. All Rights Reserved.  
Graham-Field is a trademark of GF Health Products, Inc.