

## PATIENT ROOM

# OVERBED TABLE

The Lumex® Overbed Table combines an attractive, ergonomic design with features that reduce infection control issues by minimizing areas in which dust and dirt can accumulate Manufactured with the lowest base in the industry, our Overbed tables move easily under most hospital beds or recliners.

#### **Features**:

- Lowest profile base on the market
- Exposed Graphite Grey steel U base with no plastic cover allows for easy cleaning of base and casters
- Laminate or Thermofoil top that provides durability
- Easy, one-hand activated secondary surface reduces infection control issues and eliminates hidden spaces between surfaces for bacteria to collect
- Integrated accessories
- Integrated cup holders, and spill edge<sup>2</sup>
- Optional Vanity drawer with mirror slides out from either side of the table
- Shock-absorbent gas spring cylinder for easy height adjustment
- Two (2) Year Limited Warranty



### Specifications:

	<sup>Overall</sup> Width		<sup>Overall</sup> Depth		Overall Hei low		hight Range		Main Top Width		Flip Top Width		Top Depth		Base Width		Base Depth		Base Profile Height		Product <sup>®</sup> Weight		Weight <sup>1</sup> Capacity evenly distributed	
	inch	cm <sup>+</sup>	inch	cm <sup>‡</sup>	inch	cm <sup>+</sup>	inch	cm <sup>+</sup>	inch	cm <sup>†</sup>	inch	cm <sup>+</sup>	inch	cm <sup>+</sup>	inch	cm <sup>‡</sup>	inch	cm <sup>+</sup>	inch	cm+	lb	kg‡	lb	kg⁺
Single Top	32.50	82.55	19.00	48.26	28.00	71.12	44.00	111.76	32.00	81.28	—	—	17.00	43.18	32.50	82.55	19.00	48.26	1.94	4.93	34	15.42	100	45.36
Flip Top	46.00	116.84	19.00	48.26	28.00	71.12	44.00	111.76	32.00	81.28	13.50	34.29	17.00	43.18	32.50	82.55	19.00	48.26	1.94	4.93	41	18.60	100	45.36

+ standard rounding method used 1 maximum weight capacity - evenly distributed - refer to user manual 2 thermofiol tops only 8 approximate measurement, options may alter spec

GF Health Products, Inc. is an intertek ISO 13485:2016 / MDSAP Certified Company

PColor variations may occur due to monitor or printer settings.

#### +1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. www.grahamfield.com



© 2019, GF Health Products, Inc. All Rights Reserved. Lumex is a registered trademark of GF Health Products, Inc

