

FR565WG

LUMEX[®]
by graham-field

The Extra Wide Preferred Care Recliner is ideal for clinical and acute care settings. This recliner allows for almost infinite recline positioning. The unique position lock can be caregiver activated to stabilize a position for patient comfort. A unique suspension seating system self-adjusts to each person, cradling the resident and virtually eliminating the risk of "bottoming out." The quick release upholstery allows for easy cleaning, removal and replacement. The reduced armrest and seat height enables the recliner to fit under dining and activity tables.



EXTENDED
CARE



CLINICAL
CARE



**400
POUND**¹
WEIGHT CAPACITY
EVENLY DISTRIBUTED

PACU	Rehab	Oncology	IV Therapy	Exam Room	Patient Room	Emergency Room	Dialysis (Facility)	Dialysis (In Home)	Nocturnal Dialysis	Physical Therapy	Pain Management	Post Operative Care	Same Day Surgery / ASC
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
AREAS OF USE													

HEALTHCARE SEATING SOLUTIONS

graham
field

Model

FR565WGxxx

† xxx REPRESENTS FABRIC COLOR

 **FR565WG** Product Info



Positions: INFINITE POSITIONING + CAREGIVER -ACTIVATED POSITION LOCK

Infinite Recline Positions: Allows Patient to Set a Position That is Most Comfortable
Recline Position Lock: Allows Caregiver to Secure the Chair in Any Desired Position



Features:

- Unique Ultra-flex system cradles resident to greatly reduce the risk of "bottoming out"
- 24 inch wide seat provides excellent patient comfort
- Quick-release upholstery allows for easy cleaning, removal and replacements
- Wide and soft armrest for patient comfort
- Vinyl upholstery is sulfide stain, mildew and oil resistant
- 5 inch Tente® swivel casters for ease of maneuverability: 3 Total Locking and 1 Directional Locking
- Durable epoxy coated steel frame is easy to clean and maintain
- Wall Saver armrests increase durability and keep walls from marring

Accessories:

★ sold separately

item no.†

Headrest Pillow	FR5659xxxUS
Adjustable Footrest	LX5640G
IV Pole	2101A
IV Mount	5782RG



RECLINE POSITION LOCK



SOFT, WIDE ARMREST



TENTE® CASTERS



QUICK-RELEASE
UPHOLSTERY



ANTI-TIPPERS

Specifications:

Overall Width		Overall Depth		Overall Depth fully reclined		Overall Height		Arm Width		Width between arms		Arm Height from floor		Seat Width		Seat Depth		Seat Height		Back Width		Back Height above seat		Caster Diameter		Back Angle		Weight ¹ Capacity EVENLY DISTRIBUTED	
inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	°	varies	lb	kg†
32.00	81.28	37.00	93.98	66.50	168.91	48.00	121.92	4.50	11.43	24.00	60.96	28.00	71.39	24.00	60.96	19.50	49.53	20.00	50.80	24.00	60.96	30.50	77.47	5.00	12.70	110	varies	400	181.44

* Ships Freight: request additional shipping quote

† xxxxx in the product number is the placeholder for fabric choice

† standard rounding method used

¹ maximum weight capacity - evenly distributed - refer to user manual

upright
legrest extended
fully reclined



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Available Colors:

COLOR OPTIONS AND FABRIC SWATCHES ON WEBPAGE

 **HC Fabric Choices**

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.
Tente is a registered trademark of TENTE International GmbH.

**graham
field**

GF09000.07RevF25