

FR601

LUMEX®
by graham-field

The Pivot-Arm Clinical Care Recliner is our flagship recliner that focuses on being the most aesthetically pleasing, comfortable, and durable recliner for use across a variety of clinical and acute care settings. This recliner features pivoting arms for easy and convenient cleaning, provides exceptional patient comfort with the use of Viscoelastic Memory Foam in the seat and backrest, has three comfort positions: upright, first recline and second recline. The red color-coded side pedals quickly put the patient into Trendelenburg (or shock) position while allowing the caregiver to maintain eye contact with patient.



CLINICAL
CARE



**400
POUND™**
WEIGHT CAPACITY
EVENLY DISTRIBUTED

PACU	Rehab	Oncology	IV Therapy	Exam Room	Patient Room	Emergency Room	Dialysis (Facility)	Dialysis (In Home)	Nocturnal Dialysis	Physical Therapy	Pain Management	Post Operative Care	Same Day Surgery / ASC
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

AREAS OF USE

HEALTHCARE SEATING SOLUTIONS

**graham
field™**

Options:

* Model	Urethane Arm Caps	Fully Upholstered Arms	Heat and Massage	Heat Only	Extended Legrest ³
FR601Uxxx	■				
FR601UHxxx	■	■			
FR601UHOxxx	■		■		
FR601ULxxx	■				■
FR601ULHxxx	■	■			■
FR601ULHOxxx	■		■		■
FR601Pxxx		■			
FR601PHxxx		■	■		
FR601PHOxxx		■		■	
FR601PLxxx		■			■
FR601PLHxxx		■	■		■
FR601PLHOxxx		■		■	■

† xxx REPRESENTS FABRIC COLOR

FR601 Product Info

Positions: THREE POSITION + CAREGIVER-ACTIVATED TRENDLENBURG

Three Positions for Patient Comfort: Upright, First Recline and Second Recline.

Fourth Position: Trendelenburg (or Shock) Position, Caregiver-Activated Using Side Foot Pedal on Either Side



Features:

- Both arms pivot outward, enabling easier cleaning and access
- Aesthetically pleasing design
- Wide contoured seat and back provides an enhanced level of patient comfort
- Viscoelastic Memory Foam engineered into the seat and backrest for ultimate comfort
- 5 inch Tente® swivel casters for ease of maneuverability: 3 Total Locking and 1 Directional Locking
- Two self-storing, fold down, side tables with smooth surface for easy cleaning
- Urethane armrest caps are comfortable for patient and easy to maintain; also available with fully upholstered arms
- Pullout footplate³ for easier patient entry and exit and safer patient transport
- Rear push handle for improved patient transport




FOLD DOWN, SIDE TABLES



PULLOUT FOOTPLATE³



TENTE® CASTERS



CAREGIVER ACTIVATED FOOT PEDALS



REAR PUSH HANDLE

Specifications:

Overall Width		Overall Width side tables up		Overall Depth		Overall Depth trendelenburg		Overall Height		Arm Width		Width between arms		Arm Height from floor		Seat Width		Seat Depth		Seat Height		Back Width		Back Height above seat		Back Angle			Weight ¹ Capacity	
inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	upright	legrest extended	fully reclined	lb	kg†
34.50	87.63	54.00	137.16	38.00	96.52	72.00	182.88	51.50	130.81	3.00	7.62	25.00	63.50	28.50	72.39	22.50	57.15	19.00	48.26	22.00	68.58	30.50	77.47	30.00	76.20	120	135	160	400	181.44

*Ships Freight: request additional shipping quote

†-xxx in the product number is the placeholder for fabric choice

† standard rounding method used

¹ maximum weight capacity - evenly distributed - refer to user manual

³ pullout footplate is not compatible available with extended legrest option. alters weight specs

⁵ requires factory installation include item with purchase

★ options available at additional cost - accessories sold separately

Available Colors:

COLOR OPTIONS AND FABRIC SWATCHES ON WEBPAGE HC Fabric Choices



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website

www.grahamfield.com



© 2019, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.
Tente is a registered trademark of TENTE International GmbH.

grahamfield

GF-1900031RevF25