

FR565

LUMEX[®]
by graham-field

FR565G Preferred Care Recliner allows for almost infinite recline positioning. The unique position lock can be caregiver activated to stabilize a position for patient comfort. A unique suspension seating system self-adjusts to each person.

FR565DG Drop-Arm Preferred Care Recliner has the same great features as FR565G with an easy-to-actuate self-storing drop-arms for a safer and easier lateral transfer for the resident and caregiver.

FR565TG Tilt-in-Space Drop-Arm Preferred Care Recliner features self-storing, height-adjustable drop-arms for safer, easier lateral transfer. The Tilt-in-Space function allows the caregiver to "Tilt" the seat and back for pressure relief and position change.



Drop-Arm model
FR565DG shown



Tilt-in-Space model
FR565TG shown

EXTENDED
CARE



HOME
CARE



250
POUND[†]
WEIGHT CAPACITY
EVENLY DISTRIBUTED

AREAS OF USE	PACU	Rehab	Oncology	IV Therapy	Exam Room	Patient Room	Emergency Room	Dialysis (Facility)	Dialysis (In Home)	Nocturnal Dialysis	Physical Therapy	Pain Management	Post Operative Care	Same Day Surgery / ASC
	■						■				■			

HEALTHCARE SEATING SOLUTIONS

graham
field

Options:

* Model	Recliner	Drop-Arm	Tilt-in-Space
FR565Gxxx	■		
FR565DGxxx	■	■	
FR565GTGxxx	■	■	■

† xxx REPRESENTS FABRIC COLOR

[FR565G Product Info](#)

[FR565DG Product Info](#)

[FR565TG Product Info](#)



FR565G



FR565DG



FR565TG

Positions: INFINITE POSITIONING + CAREGIVER-ACTIVATED POSITION LOCK

Infinite Recline Positions: Allows Patient to Set a Position That is Most Comfortable
 Recline Position Lock: Allows Caregiver to Secure the Chair in Any Desired Position

FR565TG: + CAREGIVER-ACTIVATED TILT-IN-SPACE POSITION

Tilt-In-Space: Caregiver-Activated Side Rocker Pedal "Tilts" Seat and Back for Pressure Relief and Position Change



Features:

- Tilt-in-Space function allows caregiver to tilt the seat and back for pressure relief and position change (TG)
- Innovative, self-storing, easy to actuate drop-arms for safer lateral transfer (DG & TG)
- Unique Ultra-flex system cradles resident to greatly reduce the risk of "bottoming out"
- Quick-release upholstery allows for easy cleaning, removal and replacements
- Vinyl upholstery is sulfide stain, mildew and oil resistant
- Lower armrest and seat height to fit under dining/activity tables
- Wide and soft armrest for patient comfort
- 5 inch casters for easy transport: 2 Rear Swivel and Wheel Locking; 2 Fixed Front Wheels
- Durable epoxy coated steel frame is easy to clean and maintain
- Wallsaver armrests increase durability and keep walls from marring

Accessories:

* sold separately	item no.†
Headrest Pillow	FR5659xxxJS
Head Support	FR56461xxx
Body Bolster	FR56471xxx
Adjustable Footrests	LX5640G
Foot-Drop Assembly	5657Gxxx
Activity Tray Table	5644G
Anti-Tip Wheels	LX5656G
IV Pole	2101A
IV Mount	5782RG



EPOXY COATED STEEL FRAME



SOFT, WIDE ARMREST



5 INCH CASTERS



QUICK-RELEASE UPHOLSTERY



REMOVABLE SIDE PANELS




SELF-STORING DROP-ARMS

Specifications:


Series	Overall Width		Overall Depth		Overall Depth fully reclined		Overall Height		Arm Width		Width between arms		Arm Height drop-arm up from floor		Seat Width		Seat Depth		Seat Height		Back Width		Back Height above seat		Back Angle	Weight ¹ Capacity			
	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†		upright	legrest extended	fully reclined	lb
FR565G					66.00	167.64					19.00	48.26	26.00	66.04					20.00	50.80					110	varies	150	250	113.40
FR565DG	28.00	71.12	36.00	91.44			48.00	121.43	4.50	11.43	19.50	49.53	26.75	67.95	19.00	48.26	19.50	49.53	22.00	55.88	19.00	48.26	31.00	78.74					
FR565TG					—	—																		110	—	—			

* Ships Freight: request additional shipping quote † xxx in the product number is the placeholder for fabric choice † standard rounding method used
 † maximum weight capacity - evenly distributed - refer to user manual * available at additional cost - accessories sold separately

Available Colors:



427 BLUE RIDGE



857 JADE



863 ROSEWOOD



GF Health Products, Inc. is an
 intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
 The most current and complete product information can be found on our website.
www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
 Lumex is a registered trademark of GF Health Products, Inc.
 Tente is a registered trademark of TENTE International GmbH.

graham
 field

GF-2400161RevA25