

FR587W

LUMEX[®]
by graham-field

The Extra Wide Clinical Care Recliner is ideal for clinical, acute care, and bariatric settings. This recliner provides patients more space than the other recliners while providing excellent value that maximizes patient comfort, product durability, and ease of use. It features patient comfort and support in all three positions: upright, first recline and second recline. Plus Trendelenburg position activated by a rear mounted release lever with dual gas cylinders for smooth assist during movement. Fold down side tables with smooth surfaces and removable side panels for easy cleaning.



BARIATRIC CARE



CLINICAL CARE



450 POUND[†]
WEIGHT CAPACITY
EVENLY DISTRIBUTED

PACU	Rehab	Oncology	IV Therapy	Exam Room	Patient Room	Emergency Room	Dialysis (Facility)	Dialysis (In Home)	Nocturnal Dialysis	Physical Therapy	Pain Management	Post Operative Care	Same Day Surgery / ASC
■	■	■	■	■	■	■	■	■	■	■	■	■	■

AREAS OF USE

Options:

* Model	Wallsaver Armrest Caps	Fully Upholstered Arms	Heat and Massage	Heat Only
FR587Wxxx	■			
FR587WHxxx	■		■	
FR587WHOxxx	■			■
FR587WxxxPA		■	■	
FR587WHxxxPA		■	■	
FR587WHOxxxPA		■		■

‡ xxx REPRESENTS FABRIC COLOR

FR587W Product Info



WALLSAVER ARM CAPS



FULLY UPHOLSTERED ARMS



Positions: THREE POSITION + CAREGIVER-ACTIVATED TRENDELENBURG

Three Positions for Patient Comfort: Upright, First Recline and Second Recline
Fourth Position: Trendelenburg (or Shock) Position, Caregiver-Activated Using Rear Mounted Hand Trigger



Features:

- Extra-Wide seat and back provides an enhanced level of patient comfort
- Higher seat & back height allows caregiver to more comfortably assist patient
- 5 inch Tente® swivel casters for ease of maneuverability: 3 Total Locking and 1 Directional Locking
- Wide and soft armrest for patient comfort
- Wallsaver armrest caps increase durability and protect walls from marring; or available with fully upholstered arms
- Two self-storing, fold down, side tables with smooth surface for easy cleaning
- Dual gas cylinders give smooth assist for Trendelenburg movement
- Smooth and removable side panels for easy cleaning and maintenance
- Rear push handle for improved patient transport



FOLD DOWN, FLUSH-MOUNT SIDE TABLES



SOFT, WIDE ARMREST



TENTE® CASTERS



REMOVABLE SIDE PANELS



REAR PUSH HANDLE

Accessories:

★ sold separately	item no.†
Headrest Pillow	FR5659xxxUS
Headrest Skirt	HRSxxx
Headrest Cover	HRC587xxx
4 Inch Extension Legs ⁶	FR5699G
Swing-away Folding Table	SAFT10
IV Pole	2101A
IV Mount	5782RG

Specifications:

Overall Width		Overall Width side tables up		Overall Depth		Overall Depth trendelenburg		Overall Height		Arm Width		Width between arms		Arm Height from floor		Seat Width		Seat Depth		Seat Height		Back Width		Back Height above seat		Back Angle upright legrest extended fully reclined			Weight ¹ Capacity EVENLY DISTRIBUTED	
inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	120	125	135	lb	kg†
33.00	83.82	53.50	135.89	45.50	115.57	78.00	198.12	50.00	127.00	4.50	11.43	23.25	59.06	28.00	71.39	24.00	60.96	20.50	52.07	22.00	55.88	30.50	77.47	33.50	85.09				450	204.12

*Ships Freight: request additional shipping quote ‡ xxx in the product number is the placeholder for fabric choice

† standard rounding method used † maximum weight capacity - evenly distributed - refer to user manual

⁶ adjust height specs: +4 inches (+10.16 cm)

★ options available at additional cost - accessories sold separately

Available Colors:

COLOR OPTIONS AND FABRIC SWATCHES ON WEBPAGE HC Fabric Choices



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.
Tente is a registered trademark of TENTE International GmbH.

graham
field

GF090005RevF25