

FR587WD

LUMEX®
by graham-field

The Extra Wide Drop-Arm Clinical Care Recliner is ideal for clinical, acute care, and bariatric settings. This recliner provides patients more space than the other recliners while providing excellent value that maximizes patient comfort, product durability, and ease of use. The drop-arms provide easier and safer lateral transfers of patients.



*shown with Bed Socket IV Pole
and IV Pole Attachment accessories*

**BARIATRIC
CARE**



**CLINICAL
CARE**



**450
POUND¹**
WEIGHT CAPACITY
EVENLY DISTRIBUTED

PACU	Rehab	Oncology	IV Therapy	Exam Room	Patient Room	Emergency Room	Dialysis (Facility)	Dialysis (In Home)	Nocturnal Dialysis	Physical Therapy	Pain Management	Post Operative Care	Same Day Surgery / ASC
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
AREAS OF USE													

HEALTHCARE SEATING SOLUTIONS

**graham
field™**

Options:

*Model	Wallsaver Armrest Caps	Fully Upholstered Arms	Heat and Massage	Heat Only
FR587WDxxx	■			
FR587WDHxxx	■		■	
FR587WDH0xxx	■			■
FR587WDxxxPA		■		
FR587WDHxxxPA		■	■	
FR587WDH0xxxPA		■		■

† xxx REPRESENTS FABRIC COLOR

[FR587WD Product Info](#)



WALLSAVER ARM CAPS



FULLY UPHOLSTERED ARMS



Positions: THREE POSITION + CAREGIVER-ACTIVATED TRENDLENBURG

Three Positions for Patient Comfort: Upright, First Recline and Second Recline

Fourth Position: Trendelenburg (or Shock) Position, Caregiver-Activated Using Rear Mounted Hand Trigger



Features:

- Innovative, self-storing, easy to actuate drop-arms for safer lateral transfer
- Extra-Wide seat and back provides an enhanced level of patient comfort
- Higher seat & back height allows caregiver to more comfortably assist patient
- 5 inch Tente® swivel casters for ease of maneuverability: 3 Total Locking and 1 Directional Locking
- Wide and soft armrest for patient comfort
- Wallsaver armrest caps increase durability and protect walls from marring; or available with fully upholstered arms
- Two self-storing, fold down, side tables with smooth surface for easy cleaning
- Dual gas cylinders give smooth assist for Trendelenburg movement
- Removable upholstered side panels for easy cleaning and maintenance
- Rear push handle for improved patient transport



FOLD DOWN, FLUSH-MOUNT
SIDE TABLES



DROP-ARM RELEASE



TENTE® CASTERS



UPHOLSTERED REMOVABLE
SIDE PANELS



REAR PUSH HANDLE

Accessories:

★ sold separately	item no.†
Headrest Pillow	FR5659xxxUS
Headrest Skirt	HRSxxx
Headrest Cover	HRC587xxx
Swing-away Folding Table	SAFT10
IV Pole	2101A
IV Mount	5782RG

Specifications:

Overall Width		Overall Width side tables up		Overall Depth		Overall Depth Trendelenburg		Overall Height		Arm Width		Width between arms		Arm Height drop-arm up from floor		Seat Width		Seat Depth		Seat Height		Back Width		Back Height above seat		Back Angle			Weight ¹ Capacity EVENLY DISTRIBUTED	
inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	upright	leg extended	fully reclined	lb	kg†
35.75	90.81	58.25	147.96	44.00	111.76	78.00	198.12	50.00	127.00	4.50	11.43	26.50	67.31	27.50	69.85	24.00	60.96	20.50	52.07	22.00	55.88	30.50	77.47	33.50	85.09	120	125	135	450	204.12

*Ships Freight: request additional shipping quote

† xxx in the product number is the placeholder for fabric choice

† standard rounding method used

¹ maximum weight capacity - evenly distributed - refer to user manual

★ options available at additional cost - accessories sold separately

Available Colors:

COLOR OPTIONS AND FABRIC SWATCHES ON WEBPAGE [HC Fabric Choices](#)



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

☎ 1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website

www.grahamfield.com



© 2018, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.
Tente is a registered trademark of TENTE International GmbH.

graham
field

GF-180001Rev D25