

FR588W

LUMEX®
by graham-field

Powered Bariatric Clinical Care Recliner supports up to 700 lb. A four-function hand pendant controls adjustable infinite positions for patient comfort. Legrest and backrest can be raised and lowered independently. Equipped with a rechargeable battery pack in case of a power failure. Quick-release trigger for Trendelenburg Position. Armrests drop for easier lateral transfers. Fold down side tables with smooth surfaces for easy cleaning. This recliner is ideal for the safety of care staff assisting bariatric patients along with the patients themselves.



BARIATRIC
CARE



EXTENDED
CARE



CLINICAL
CARE



PACU	Rehab	Oncology	IV Therapy	Exam Room	Patient Room	Emergency Room	Dialysis (Facility)	Dialysis (In Home)	Nocturnal Dialysis	Physical Therapy	Pain Management	Post Operative Care	Same Day Surgery / ASC
■	■	■	■	■	■	■	■	■	■	■	■	■	■
AREAS OF USE													

HEALTHCARE SEATING SOLUTIONS

graham
field

Model

FR588Wxxx

‡ xxx REPRESENTS FABRIC COLOR

FR588W Product Info

Fully Upholstered Arms



HAND PENDANT
ADJUSTABLE POSITIONS
FOR PATIENT COMFORT



shown with
Swing Away Folding Table
accessory

Positions: INFINITE POSITIONING + CAREGIVER-ACTIVATED TRENDLENBURG

Infinite Recline Positions: Allows Patient to Set a Position That is Most Comfortable

Fourth Position: Trendelenburg (or Shock) Position, Caregiver-Activated Using Rear Mounted Hand Trigger



Features:

- Fully Upholstered Arms are comfortable
- Powered Recline with four-function controller
- Legrest can be raised and lowered independently of backrest
- Battery backup in case of a power failure
- Rear push handle for improved patient transport
- Innovative, self-storing, easy to actuate drop-arms for safer lateral transfer
- Two self-storing, fold down, side tables with smooth surface for easy cleaning
- 5 inch casters for easy transport: 2 Rear Swivel and Wheel Locking; 2 Fixed Front Wheels

Accessories:

★ sold separately

Swing Away Folding Table

item no.†

SAFT30



POWERED LEGREST



EASY TO ACTUATE
DROP-ARMS



5 INCH CASTERS



FOLD DOWN, FLUSH-MOUNT
SIDE TABLES



REAR PUSH HANDLE



HEADREST PILLOW
INCLUDED

Specifications:

Overall Width		Overall Width <i>side tables up</i>		Overall Depth		Overall Depth <i>trendelenburg</i>		Overall Height		Arm Width		Width <i>between arms</i>		Arm Height <i>from floor</i>		Seat Width		Seat Depth		Seat Height		Back Width		Back Height <i>above seat</i>		Back Angle	Weight ¹ Capacity <small>EVENLY DISTRIBUTED</small>		
inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†		lb	kg†	
38.00	96.52	48.00	121.92	59.00	149.86	77.00	195.58	53.00	134.62	5.00	12.70	24.50	62.23	31.00	78.74	26.50	67.31	23.50	59.69	25.50	64.77	32.00	81.28	35.00	88.90	110	varies	700	317.52

* Ships Freight: request additional shipping quote

‡ xxx in the product number is the placeholder for fabric choice

† standard rounding method used

¹ maximum weight capacity - evenly distributed - refer to user manual

★ accessories sold separately

Available Colors:



427 BLUE RIDGE

851 DOE SKIN

857 JADE

863 ROSEWOOD



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.
Tente is a registered trademark of TENTE International GmbH.

graham
field

GF0900052RevF25