


WHEELCHAIR ACCESSORIES

CUP HOLDER

Features

- Easily snaps and adjusts by tension clip-ons onto walkers, rollators, wheelchairs, transport chairs, walkers, quad and standard canes with tubing diameter from 5/8" to 1" (1.58 cm to 2.54 cm)
- Tension adjustable clamp bottle holder united to a 3 1/2" (8.89 cm) diameter shelf base by two screws securing two metal clip-ons for extra support
- Frame made with durable high-density polypropylene
- 3 1/2" (8.89 cm) diameter composite shelf base to support bottle
- Total product weight: 2.6 ounces (73.70 g)
- Product dimensions: 3" x 3 1/2" x 3" (7.62 cm x 8.89 cm x 7.62 cm)
- Fits 2 1/2"-3" (6.35 cm x 7.62 cm) in diameter bottles (holder can be adjusted to fit bottle size)
- Weight capacity: bottled beverage containers from 8 to 20 oz (236.59 ml to 591.47 ml) evenly distributed – Maximum diameter 2.5" to 3" (6.35 cm x 7.62 cm)
- Warranty: 6 (six) months limited warranty.



	overall Width		overall Depth		overall Height		maximum Diameter		Product Weight		UOM
	inch	mm	inch	mm	inch	mm	inch	mm	oz	g	
EJCUPH19	3	7.6	3.5	8.9	3	7.6	2.5-3	6.35-7.6	2.6	73.7	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved.
Everest & Jennings is a registered trademark of GF Health Products, Inc.

graham
field

GF2000146RevB25