

BARIATRIC/HEAVY DUTY WHEELCHAIRS

PARAMOUNT XD WHEELCHAIR

The Paramount XD is an alternative for the bariatric user who prefers a manual wheelchair. It features 26" and 30" seat widths as well as a weight capacity of 650 lb. Meets Medicare specifications for K0007 for wheelchairs.

Features

- Available in 26" and 30" seat widths with 19 ½" seat depth
- Black reinforced nylon upholstery
- Standard chart pocket on back upholstery for models after serial number AIS338WC0046-0055
- Reinforced welds and gussets on frame with double crossbrace
- Two-position axle and caster fork provides 2" seat height adjustment
- 8"x 2" non-marking solid casters and 24" solid rear wheels with chromed handrims
- Detachable padded desk arms for easy transfer
- Adjustable swingaway footrest or elevating legrest with extra wide metallic footplate
- Maximum weight capacity 650 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0007

Warranty:

Five (5) year limited warranty



desk arms (shown) and
swingaway footrest (shown)
— other configurations available



with Elevating Legrests















Chart pocket



Color variations may occur due to monitor or printer settings.



Paramount XD Heavy Duty Wheelchair

Item	HCPCS [‡] code	Overall Chair Width		Overall Chair Depth		Overall Chair Height		Folded Chair Width		Seat Width		Seat Depth		Seat Height		Back Height		Wheel Diameter		Caster Diameter		Product [§] Weight w/o front rigging		Weight [†] Capacity EVENLY DISTRIBUTED		Detach. Full Detach. Desk Swingaway Footrest Elevating Legrest Rear Anti-Tippers							
																																	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg [†]	lb	kg [†]								
5PX10620	K0007/E2202	35.5	90.2	49	124.5	39	99.1	14.5	36.8	26	66	19.5	49.5	18-20	45.7-50.8	18.5	47	24	61	8	20.3	93	42.2	650	294.8	—	■	■	—	■	■		
5PX10630																																	
5PX10820	K0007/E2202	40	101.6	49	124.5	39	99.1	14.5	36.8	30	76.2	19.5	49.5	18-20	45.7-50.8	18.5	47	24	61	8	20.3	93	42.2	650	294.8	—	■	■	■	—	■	■	
5PX10830																											—	■	—	■	—	■	■

■ Standard

■ Optional

■ Available

■ Accessory

— Not Apply / Available

‡ designed to meet the requirements of

† standard rounding method used

† maximum weight capacity - evenly distributed - refer to user manual

[‡] designed to meet the requirements of
 [†] standard rounding method used
 [†] maximum weight capacity - evenly distributed - refer to user manual
[§] approximate measurement: product weight — as seat width and rigging change, so does product weight
 ★ options available at additional cost - accessories sold separately

Parts & Accessories: ★ sold separately

	HCPCS [‡]	UOM	Item No.
Wheel lock extension, for Paramount XD	E0961	pair	90762936
Positioning belt for Paramount XD (metal buckle)	E0978	each	90762937
Anti-folding bar, chromed (lower frame) with telescopic width (22" - 30")	—	each	92523851

	HCPCS [‡]	UOM	Item No.
Anti-tippers, for Paramount XD	—	pair	9A000022
Replacement Knobs	—	pair	90762938
Combination O2 & IV Holder	E2208	each	B10320F

	HCPCS [‡]	UOM	Item No.
O2 Cylinder Holder	E2208	each	B10325
Wheelchair Positioning Tray	—	each	EJTRAY18



Anti-folding bar



Positioning belt



Anti-tippers



GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2019, GF Health Products, Inc. All Rights Reserved.
Everest & Jennings is a registered trademark of GF Health Products, Inc.

grahamfield

GF0700185RevH25