

RECLINING WHEELCHAIRS

ADVANTAGE® RECLINER

The Advantage Recliner comes with a higher back and removable full support headrest. The headrest includes an additional neck support that adjusts easily to fit every need. The durable hammertone frame is designed to support up to 300 lb and in the bariatric version up to 450 lb users.

Features


- Available in 16", 18" and 20" seat widths, with 17" seat depth and in 22" seat width with 19" seat depth
- Features manual full reclining back range from 90-170 degrees
- Removable full support headrest
- Detachable padded desk or full arms
- Hammertone finish provides durability and style
- Embossed black leatherette type vinyl upholstery for easy cleaning
- Anti-tippers standard
- Comes with spreader metallic back bar which works as anti-fold device
- Industry compatible tool-free adjustable elevating front rigging
- Standard chart pocket on back
- Maximum weight capacity 300 lb - EVENLY DISTRIBUTED for 16,18 and 20" models
- Maximum weight capacity 450 lb - EVENLY DISTRIBUTED for 22" models
- 16", 18", 20" models meet Medicare code K0001/ E1226 (Recliner Back) + E0966 (Head Rest Extension) + E0971
- (Anti-tipping device) + K0195 (Elevating Legrests)
- 22" models meet Medicare code K0007/E1226
- (Recliner Back + E0966 (Head Rest Extension) + E0971
- (Anti-tipping device) + K0195 (Elevating Legrests) + E2201

Warranty:

Five (5) year limited warranty



desk arms (shown) and
elevating legrest (shown)
— other configurations available





Advantage® Recliner

Item	HCPCS [‡] code	Overall Chair Width		Overall Chair Depth		Overall Chair Height		Folded Chair Width		Seat Width		Seat Depth		Seat Height		Back Height		Wheel Diameter		Caster Diameter		Product [§] Weight w/o front rigging		Weight [†] Capacity EVENLY DISTRIBUTED		Detach. Full arms	Detach. Desk	Swingaway Footrest front rigging	Elevating Legrest	Rear Anti-Tippers
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg†	lb	kg†					
3K010130	K0001/ E1226	26.5	67.3	34	86.4	51	129.5	11.5	29.2	18	45.7	17	43.2	19.5	49.5	23	58.4	24	61	8	20.3	53	24	300	136.1	—	■	—	■	■
3K010150																										■	—	—	■	■
3K010230	K0001/ E1226	28.5	72.4	34	86.4	51	129.5	11.5	29.2	16	40.6	17	43.2	19.5	49.5	23	58.4	24	61	8	20.3	53	24	300	136.1	—	■	—	■	■
3K010250																										■	—	—	■	■
3K010330	K0001/ E1226	30.5	77.5	34	86.4	51	129.5	11.5	29.2	20	50.8	17	43.2	19.5	49.5	23	58.4	24	61	8	20.3	53	24	300	136.1	—	■	—	■	■
3K010350																										■	—	—	■	■
3K020430	K0007/E1226	32.5	82.6	34	86.4	51	129.5	11.5	29.2	22	55.9	19	48.3	19.5	49.5	23	58.4	24	61	8	20.3	57	25.9	450	204.1	—	■	—	■	■
3K020450																										■	—	—	■	■

■ Standard

■ Optional

■ Available

■ Accessory

— Not Apply / Available

† designed to meet the requirements of

† standard rounding method used

† maximum weight capacity - evenly distributed - refer to user manual

■ Standard ■ Optional ■ Available ■ Accessory — Not Apply / Available

‡ designed to meet the requirements of † standard rounding method used † maximum weight capacity - evenly distributed - refer to user manual
§ approximate measurement: product weight — as seat width and rigging change, so does product weight ★ options available at additional cost - accessories sold separately

Parts & Accessories: ★ sold separately

	HCPCS [‡]	UOM	Item No.
Footrest Heel Loop	E0951	pair	90763064
Wheel Lock Lever Extension	E0961	pair	90763062
Positioning Belt (auto-style, push button) (adjustable length: 27 - 47 inches (68.58 - 119.38 cm))	E0978	each	300-505
Anti-tippers, for 16x17/18x17/20x17 (300 lb)	E0971	pair	90763730
Anti-tippers, for 22x19 (450 lb)	E0971	pair	90763067
Anti-Fold Bar (Chromed) Lower Frame	—	each	92523850

	HCPCS [‡]	UOM	Item No.
Hammertone back spreader bar 16" (reclining back)	—	each	907637E1
Hammertone back spreader bar 18" (reclining back)	—	each	907637E2
Hammertone back spreader bar 20" (reclining back)	—	each	907637E3
Hammertone back spreader bar 22" (reclining back)	—	each	907637E4
Headrest posts,	—	pair	90763720

	HCPCS [‡]	UOM	Item No.
Headrest cushion, black vinyl	—	each	90763735
Headrest upholstery black 16" (includes hardware)	—	each	907637D1
Headrest upholstery black 18" (includes hardware)	—	each	907637D2
Headrest upholstery black 20" (includes hardware)	—	each	907637D3
Headrest upholstery black 22" (includes hardware)	—	each	907637D4



GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2019, GF Health Products, Inc. All Rights Reserved.
Everest & Jennings is a registered trademark of GF Health Products, Inc.

graham field

GF1200130RevG25