

EXAM ROOM SUPPLIES

STAINLESS STEEL GRADUATED MEDICINE CUP

Features

- Used to measure liquids
- Capacity: 2 oz, 60 cc, or 4 tbsp
- Labeled in ounces, cc (cubic centimeters) and tablespoons
- Autoclavable
- 18-8 (304) stainless steel with seamless construction



	overall Height		overall Width		Volume			UOM
	inch	cm	inch	cm	oz	tbsp	ml	
3241	2	5.1	1.75	4.4	2	4	59.1	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2009, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.