

NOW IS THE TIME TO REFRESH YOUR OFFICE FOR A NEW START IN 2025

Hausted[®] by graham-field TREATMENT TABLES

Hausted treatment tables are designed to provide a comfortable, and supportive base for treatment and diagnostics in an array of medical settings. These tables are available in static or adjustable models with or without pass through drawers.



500 POUND
WEIGHT CAPACITY,
EVENLY DISTRIBUTED

4100
SERIES
flat top

‡4100-xx shown in aqua green | AQ



4103
SERIES

shelf, 2 pass-through drawers, and flat top
‡4103-DD-xx shown in dusty blue | DB



4105
SERIES

shelf, adjustable back and contour top
‡4105-xx shown in pearl | PL



4105
SERIES
shelf, 2 pass-through drawers, adjustable back and contour top
‡4105-DD-xx†

Hausted[®] by graham-field EXAM ROOM TREATMENT TABLE

All 4200 series exam tables have three metal pass-through drawers rather than the typical plastic, so items can be accessed on both sides

- Two Storage and Three Pass Through Drawers
- Pull-Out Leg Rest with Integrated Treatment Tray
- Five Position Continuous Lock Stirrups
- All-metal base table construction
- Powder Coat Finish - Foghorn
- Pelvic Tilt

500 POUND
WEIGHT CAPACITY,
EVENLY DISTRIBUTED



4200
SERIES
4200-xx†

shown in saddle tan | ST

Options: ‡4200-E0-xx Integrated Outlet
*‡4200-DW-xx Drawer Warmer + Integrated Outlet

Hausted[®] by graham-field HIGH / LOW POWER EXAM TABLE

The medical exam table is essential in any health care facility. Hausted provides high-quality options to meet the needs of hospitals, clinics, ERs and more.

- An all-metal construction for the base table provides longevity
- Pull-Out Leg Rest with Integrated Treatment Tray
 - ◊ Tray is independent of the Legrest for any time, easy access
- Hand Pendant controls hi-low, back, flat, return to chair
- UL listed (60601-1) DC electromechanical drive
- Five Position Continuous Lock Stirrups
- Powder Coat Finish - Foghorn
- Power Assisted Back
- Pelvic Tilt

600 POUND
WEIGHT CAPACITY,
EVENLY DISTRIBUTED

4600
SERIES

‡4600-xx

shown in pearl | PL
and in artichoke | AR

Options: ‡4600-FS-xx Foot Switch
*‡4600-AA-xx Assist Arms
‡4600-PS-xx Storage Drawer
‡4600-PS-DW-xx Drawer with Warmer
‡4600-FL-xx Fully Loaded with All Options‡



Hausted[®] by graham-field REPLACEMENT TOPS

Make your exam tables look new at a fraction of the cost. Our replacement tops fit many major brand tables.



100

‡GF100RT-xx
shown in violet | VI

104

‡GF104RT-xx
shown in teal | TE

204

‡GF204RT-xx
shown in british tan | BT

304

‡GF304RT-xx
shown in ocean gray | OG

404

‡GF404RT-xx
shown in doe | DO

405

‡GF405RT-xx
shown in shadow | SH

420

‡GF420RT-xx
shown in artichoke | AR

460

‡GF460RT-xx
shown in british tan | BT



EXAM ROOM FURNISHINGS

A COMPLETE EXAM ROOM SOLUTION

Lumex® products complete the Hausted® Exam Room Furnishings line and are available in color matched vinyl upholstery. DEHP-free vinyl covering available in 24 industry standard colors. See our online swatch card here [📄](#)

LUMEX® by graham-field BARIATRIC PHLEBOTOMY CHAIRS

The Lumex 6200 Series bariatric blood draw/phlebotomy chair is designed with both the patient and phlebotomist in mind. The sensible and affordable design offers a small footprint and modern style while keeping patient comfort and well-being the highest priority.

- Fully upholstered padded seat, back, and armrests
- Powder coat frame finish - Foghorn
- Individual height adjustable arm rests
- Single Flip Arm standard
- Dual Flip Arm option available
- Shown with wall saver wireframe accessory[§] (6200-WALL SAVER)[‡]



6200
SERIES
Single Flip Arm
shown in british tan | BT
‡6200-xx



6200
SERIES
Dual Flip Arm
shown in blue | BU
‡6200DA-xx[‡]

Hausted® by graham-field SIDE CHAIRS

Looking for comfortable and contemporary patient seating for your lobby or waiting room? Our wall saver design keeps the backrest further away from the wall to reduce scratching or marring keeping maintenance costs down and your facility looking its best.

- 14 Gauge Steel Tubing
- Powder Coat Finish-Black
- Wall Saver Leg Design



2000
SERIES
Wall Saver Side Chair
shown in violet | VI
‡2000-xx



2010
SERIES
Wall Saver Arm Chair
shown in violet | VI
‡2010-xx



2000
SERIES
Wall Saver Arm Chair
Slant Arms
shown in violet | VI
‡2020-xx

LUMEX® by graham-field PHYSICIAN STOOLS

We offer a variety of medical stools to meet the needs of an active facility. Offering preferred stool features manufactured in the USA with high quality materials to ensure a long life.

9000 SERIES Physician Exam Stool

- 16 inch (40.64 cm) seat diameter
- Choice of three base styles
- Choice of three hand controlled pneumatic assist height adjustment



9501
SERIES
Black Composite Base
5 configurations available
shown in gray | GR
‡9501-xx



9702
SERIES
Toe Cap Aluminum Base
5 configurations available
shown in plum | PU
‡9702-xx



9903
SERIES
Polished Aluminum Base
3 configurations available
shown in british tan | BT
‡9903-xx



7408
SERIES
Brushed Aluminum Base
3 configurations available
shown in pearl | PL
‡7408-xx

7000 SERIES Physician Surgical Stool

- 16 inch (40.64 cm) seat diameter
- Brushed aluminum base
- Foot-controlled hydraulic height adjustment



Options:
*†

with Backrest **9501-B-xx**
with Footring **9501-FR-xx**
with Backrest & Footring **9501-BFR-xx**
Lab Height **9501-10BFR-xx**
with Backrest and Footring

with Backrest **9702-B-xx**
with Footring **9702-FR-xx**
with Backrest & Footring **9702-BFR-xx**
Lab Height **9702-10BFR-xx**
with Backrest and Footring

with Footring **9903-FR-xx**
Lab Height with Footring **9903-10FR-xx**

with Backrest **7408-B-xx**
with Crescent Backrest **7408-CB-xx**



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

‡ xx in the product number is the placeholder for vinyl/finish choice
§ representative image shows this option * options available at additional cost
1. maximum weight capacity - evenly distributed - refer to user manual 2. accessories sold separately

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.hausted.com | www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Hausted is a registered trademark of GF Health Products, Inc.

