

**The Ortho-Biotic™ II Recliner** is a premium aesthetically pleasing recliner that provides ultimate comfort for a variety of applications in extended care, clinical care, and acute care settings. Features include, infinite recline positions, independent legrest operation, and the use of Viscoelastic Memory Foam in the seat and backrest for ultimate comfort. Using an outside hand lever, the patient can control the back angle to find the position that is most comfortable to them.



EXTENDED  
CARE



CLINICAL  
CARE



**350  
POUND**<sup>1</sup>  
WEIGHT CAPACITY  
EVENLY DISTRIBUTED

PACU	Rehab	Oncology	IV Therapy	Exam Room	Patient Room	Emergency Room	Dialysis (Facility)	Dialysis (In Home)	Nocturnal Dialysis	Physical Therapy	Pain Management	Post Operative Care	Same Day Surgery / ASC
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
AREAS OF USE													

## Options:

*Model	Wood Arm Caps	Fully Upholstered Arms	Heat and Massage	Heat Only
FR597Gxxx	■			
FR597GHxxx	■		■	
FR597GH0xxx	■			■
FR597Pxxx		■		
FR597PHxxx		■	■	
FR597PH0xxx		■		■

‡ xxx REPRESENTS FABRIC COLOR

 **FR597** Product Info



## Positions: INFINITE POSITIONING + CAREGIVER-ACTIVATED TRENDENBURG

Infinite Recline Positions: Allows Patient to Set a Position That is Most Comfortable  
Trendelenburg (or Shock) Position, Caregiver-Activated Using Side Foot Pedal



## Features:

- Aesthetically pleasing design
- Viscoelastic Memory Foam engineered into the seat and backrest for ultimate comfort
- Legrest-Recline lever activates independent movement from infinite backrest position
- Headrest pillow included for extra patient comfort
- 4 inch Tente® swivel casters for ease of maneuverability: 3 Total Locking and 1 Directional Locking
- Fully Upholstered armrest are comfortable for patient and easy to maintain; also available with wood arm caps
- Pullout footplate for easier patient entry and exit and safer patient transport
- Rear push handle for improved patient transport



HEADREST PILLOW



PULLOUT FOOTPLATE



TENTE® CASTERS



LEGREST-RECLINE LEVER



REAR PUSH HANDLE

## Accessories:

★ sold separately

	item no.†
Headrest Pillow	included
Headrest Skirt	HRSxxx
Left Single Side Table <sup>4,5</sup>	FR597GTBLMNT
Swing-away Folding Table	SAFT20
IV Pole	2101A
IV Mount <sup>5</sup>	FR597GIVMNT

## Specifications:

Overall Width		Overall Width <i>side tables up</i>		Overall Depth		Overall Depth <i>trendelenburg</i>		Overall Height		Arm Width		Width <i>between arms</i>		Arm Height <i>from floor</i>		Seat Width		Seat Depth		Seat Height		Back Width		Back Height <i>above seat</i>		<i>upright</i> <i>legrest</i> <i>fully recl</i>			Weight <sup>1</sup> Capacity <small>EVENLY DISTRIBUTED</small>	
inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	inch	cm†	Back Angle			lb	kg†
33.00	83.82	45.00	114.30	37.00	93.98	73.00	185.42	46.00	116.84	3.50	8.89	24.00	60.96	26.75	67.95	23.00	58.42	20.00	50.80	21.00	53.34	23.00	58.42	30.00	76.20	120	139	180	350	158.76

\*Ships Freight: request additional shipping quote

‡-xxx in the product number is the placeholder for fabric choice

† standard rounding method used

<sup>1</sup> maximum weight capacity - evenly distributed - refer to user manual

<sup>4</sup> supports a single table left side only

<sup>5</sup> requires factory installation include item with purchase

★ options available at additional cost - accessories sold separately

## Available Colors:

COLOR OPTIONS AND FABRIC SWATCHES ON WEBPAGE  **HC Fabric Choices**



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website

[www.grahamfield.com](http://www.grahamfield.com)



© 2014, GF Health Products, Inc. All Rights Reserved.  
Lumex is a registered trademark of GF Health Products, Inc.  
Tente is a registered trademark of TENTE International GmbH.

**graham field**

GF1400070Rev D25